

## 脂溶性ビタミンの分析

### Oil-soluble Vitamines

#### Analytical Conditions

Column	: L-column ODS 4.6mm × 150mm
Mobile phase	: Methanol/Acetonitrile =40/60
Flow rate	: 1mL/min
Temp.	: 40
Detector	: UV285nm 0.16AUFS
Sample	: Vitamin K <sub>3</sub> Vitamin A Vitamin A acetate Vitamin D <sub>2</sub> Vitamin D <sub>3</sub> Vitamin E Vitamin E acetate Vitamin K <sub>1</sub>

