

脂肪族有機酸の分析 Aliphatic organic acids

Analytical Conditions

Column	: L-column ODS 4.6mm × 250mm
Mobile phase	: 50mM NH ₄ H ₂ PO ₄ (pH=2.0)
Flow rate	: 1mL/min
Temp.	: 40
Detector	: UV210nm 0.16AUFS
Sample	: Oxalic acid Tartaric acid Glycolic acid Malic acid Malonic acid Lactic acid Maleic acid Citric acid Fumaric acid Acrylic acid

