

## 芳香族有機酸の分析

## Aromatic organic acids

### Analytical Conditions

Column	: L-column ODS 4.6mm × 250mm
Mobile phase	: 20mM H <sub>3</sub> PO <sub>4</sub> /Acetonitrile =80/20
Flow rate	: 1mL/min
Temp.	: 40
Detector	: UV254nm 0.08AUFS
Sample	: Gallic acid Protocatechuic acid o-Phthalic acid Gentisic acid m-Hydroxybenzoic acid Salicylic acid

