

ニトロフェノール類の分析

Nitrophenols

Analytical Conditions

Column	: L-column ODS 4.6mm × 150mm
Mobile phase	: 20mM NaH ₂ PO ₄ (pH=2.0)/Acetonitrile =70/30
Flow rate	: 1mL/min
Temp.	: 40
Detector	: UV254nm 0.16AUFS
Sample	: Phenol p-Nitrophenol m-Nitrophenol 2,6-Dinitrophenol 2,4-Dinitrophenol o-Nitrophenol

