

クロロフェノール類の分析

Chlorophenols

Analytical Conditions

Column	: L-column ODS 4.6mm × 250mm
Mobile phase	: 20mM NaH ₂ PO ₄ (pH=2.5)/Acetonitrile =50/50
Flow rate	: 1mL/min
Temp.	: 40
Detector	: UV254nm 0.04AUFS
Sample	: Phenol o-Chlorophenol p-Chlorophenol m-Chlorophenol 2,6-Dichlorophenol 2,4,5-Trichlorophenol 2,4,6-Trichlorophenol 2,3,4,6-Tetrachlorophenol

