


テオフィリン Theophylline


Key words : テオフィリン 医薬品
Column : USP category: L1

[Analytical conditions]

Column : L-column ODS (C18, 5 μ m, 12 nm), 4.6 mm I.D. $\times$ 150 mm L.; Cat. No. 622070
Eluent : CH₃CN/20 mM Phosphate buffer pH 6.8 (10/90)
Flow rate : 1 mL/min
Temperature : 40°C
Detection : UV 270 nm
Injection volume : -
System : -



Sample:
1. Theophylline
テオフィリン



2. Caffeine
カフェイン

